



## Westerstrand Visual count down timer

# Time Log



Westerstrand Time Log makes it possible to understand time even if one is not able to read time in a conventional manner of clocks.

Time Log illustrates notions as WHEN and FOR HOW LONG everyday routines and activities will take place.

**Time Log** is an aid for persons with cognitive dysfunction, who find our ordinary time notions too abstract. If one does not understand the ordinary clock, it is difficult to know, WHEN something will happen. There is another necessary time notion as FOR HOW LONG I have to do this.

**Time Log** is developed in cooperation with specialists for designing technical help devices. It will help to plan a day by illustrate WHEN an activity or event should start and FOR HOW LONG it should go on.

**Time Log** is available in six versions covering time periods from 2 minutes up to 8 hours. Every Tim Log's version has 4 installed time periods to choose between. Delivery includes coloured stickers corresponding to the colours of the buttons on the Time Log. By using stickers, activities and events can be scheduled as colour marks on a timetable. The count down is illustrated by a column of red diodes. The diodes will turn off one by one until all diodes are turned off. To indicate, that the selected time has expired, the whole column of diodes blinks and a short signal will sound 5 times.

**Time Log** is a battery driven device for indoor use in dry spaces. As options are available a safety-bag and keychain.



## Function

By pressing one of start buttons (red, yellow, green or blue with indication of period) corresponding to requested time, the number of diodes will turn on and count down will start. The count down is illustrated by a column of red diodes turning off one by one until all diodes are turned off. When count down has finished, the whole column of diodes will start to blink and five short signals will be given. For turning off, one of the start buttons has to be pressed.

Level of sound can be adjusted by a thumb wheel on the left side of the cover.

The time set can be erased by pressing red and blue start button simultaneously and keeping them pressed down during two (2) seconds. Thereafter a new time can be started.

**New**

Art.nr.119628-00	Art.No.119622-00	Art.No.119621-00	Art.No.119620-00	Art.No.119623-00	Art.nr.119630-00
<b>8 min 16 diodes</b>	<b>20 min 20 diodes</b>	<b>60 min 12 diodes</b>	<b>60 min 20 diodes</b>	<b>80 min 20 diodes</b>	<b>8 hours 16 diodes</b>
2 min / 4 diodes 4 min / 8 diodes 6 min / 12 diodes 8 min / 16 diodes	5 min / 5 diodes 10 min / 10 diodes 15 min / 15 diodes 20 min / 20 diodes Art.No. for option without digits – 119622-11	15 min / 3 diodes 30 min / 6 diodes 45 min / 9 diodes 60 min / 12 diodes	15 min / 5 diodes 30 min / 10 diodes 45 min / 15 diodes 60 min / 20 diodes Art.No. for option without digits – 119620-11	20 min / 5 diodes 40 min / 10 diodes 60 min / 15 diodes 80 min / 20 diodes Art.No. for option without digits – 119623-11	2 h / 4 diodes 4 h / 8 diodes 6 h / 12 diodes 8 h / 16 diodes
<b>Interval of time: 30 sec/diode</b>	<b>Interval of time: 1 min/ diode</b>	<b>Interval of time: 5 min/ diode</b>	<b>Interval of time: 3 min/ diode</b>	<b>Interval of time: 4 min/diode</b>	<b>Interval of time: 30 min/diode</b>

## Battery

The delivery includes 2 alkaline 1,5 V batteries type R6. These batteries will last for approx. 60 hours use. It is also possible to use re-chargeable batteries of type NiMH, size R6. Re-chargeable batteries are not supplied by Westerstrand.

## Data

Dimensions (L x W x H): 111 x 52 x 27 mm.  
Weight: 115 g

## Accessories/Options

Art.No.119608-00	Art.No.119607-00	Art.No.119622-10	Time Logs with 20 diodes and with front without digits	Stickers with corresponding digits for time logs with 20 diodes
<b>Safetybag with floatpillow and hook</b>	<b>Keychain with carabiner with connection loop</b>	<b>20-min Time Log with "standing" digits</b>		

## Norms

Time Log conforms with EN 61000-6-3:2005, EN 61000-4-2 och EN 61000-4-3.

**WESTERSTRAND URFABRIK AB Sweden**

P.O. Box 133, SE-545 23 TÖREBODA  
E-mail: info@westerstrand.se

Tel. +46 506-480 00 Fax +46 506-48051  
Internet: www.westerstrand.se